

Spackman Style Uterine Cannula



Simple yet effective



Spackman Style Uterine Cannula

The Femcare Uterine Cannula is based on the "Spackman" design. It is a single use device designed for manipulation of the uterus and injection of fluid.

The "Spackman" cannula is provided with an adjustable locking mechanism for holding the Vulsellum forceps in place, these are used to 'grasp' the cervix ensuring a tight seal. Once the manipulator is in place uterine manipulation is possible.

The holes in the distal end of the manipulator allow for hydrotubation to test the patency of the fallopian tubes; either by injecting dye or contrast medium and checking its flow through the uterus and fallopian tubes during Laparoscopy or Hysterosalpingography.



The unique locking device holds the Vulsellum forceps in place



Acorn to achieve an effective seal



The Vulsellum Forceps are located into the locking device and the trigger pulled back to secure the device.

Product Code	Product Description	Pack Size
SP-DIS-100	Uterine Cannula	25



Europe

Femcare Ltd.
32 Premier Way
Romsey, Hampshire SO51 9DQ
United Kingdom
Tel (UK) +44.1794.525.100
Tel (FR) +33.9.7518.5477
Fax +44.1794.525.101
enquiries@femcare.co.uk

Australia

Femcare-Australia Pty Ltd
Unit 12, 5 Gladstone Road
Castle Hill, NSW 2154
Tel (AUS) +61.2.9045.4100
Tel (NZ) +64.9.886.6884
Fax +61.2.9894.1913
customer@femcare.com.au

Canada

Femcare Canada
6355 Kennedy Road, Unit 15
Mississauga, ON L5T 2L5
Tel +1 905.795.1102
Fax +1 905.795.8154
customer@femcare.ca

Products shown may not be available in all countries.
Contact Femcare for product availability in your area.